

Cityliner₄₀₈



CITYLINER 408 AN INSPIRING INTRODUCTORY MODEL

The new Cityliner 408 provides excellent value for money, without compromise. Equipped with a 210 watt motor it allows travelling speeds of up to 4 mph, with full lights including indicators, hazards and brake lights. All these features make the new Cityliner 408 the smart way to get around.

MEYRA
ORTOPEDIA
We move people.

MODEL NO. 3.264

Cityliner₄₀₈

- 6 km/h (4 mph)
- Range up to max. 32 km (19,9 miles)
- User weight up to 130 kg (20 stone)
- Climbing ability up to 21% (with 75 kg/12 stone user weight)
- Motor power 210 W
- Adjustable steering column
- Standard-seat, height adjustable, swivelling with sliding rails and suspension
- Flip-up arm supports
- Front shopping basket
- Aluminium wheel rims



| Dimensions/size for transport | |
|--|-----------------------------|
| Length | 122 cm (48") |
| Width | 59 cm (23") |
| Seat height, adjustable | 39 – 49 cm (15,4 – 19,3") |
| Seat depth | 38 cm (14,9") |
| Seat width, adjustable | 45 – 66 cm (17,7 – 25,9") |
| Height for transport (folded tiller, without seat) | 58 cm (22.8") |
| Chassis | |
| Automatic safety brake | ● |
| Wheels | 4 |
| Wheel diameter front wheel | 26 cm (10,4") |
| Wheel diameter rear wheel | 26 cm (10,4") |
| Ground clearance | 10 cm (3,9") |
| Obstacle height, max. | 7 cm (2,8") |
| Turning radius | 1,65 m (64,9") |
| Electronic | |
| Motor constant performance | 210 W |
| Motor peak performance | 700 W |
| Batteries | 34 Ah (5 h) 40 Ah (20 h) |
| Charger | 6 A |
| Controller | 70 A |

| Performance | |
|--|---------------------|
| Travelling range, max. (under test conditions) | 32 km (19,9 miles) |
| Technical climbing ability (with 75 kg/12 stone user's weight) | 12° (21 %) |
| Safe ascending/descending gradient | 8° (14 %) |
| Weights | |
| Unladen weight | 84,5 kg (186,3 lbs) |
| Permitted user weight | 130 kg (20 stone) |
| Permitted total weight | 240 kg (38 stone) |
| Speed | |
| Speed, max. | 6 km/h (4 mph) |
| Features / Options / Accessories | |
| Colour | silvermetallic |
| Aluminium wheel rims | ● |
| Lights | ● |
| Shopping basket, additional | ● |
| Cane holder | ● |
| Integrated bumper | ● |
| Rear-view mirror, pair | ● |

● Standard ● Option

A.M. Mobility
The Solihull Mobility Centre
767-769 Old Lode Lane
Solihull
West Midlands
B92 8JE

MEYRA
ORTOPEDIA
We move people.