



Cityliner₃₁₂

Cityliner₄₁₂



CITYLINER 312 / 412 A VERSATILE SCOOTER FOR UNLIMITED MOBILITY

As a compact class 3 scooter the new Cityliner 312 and 412 are equipped with whisper quiet but powerful motor, along with all round suspension, making short or longer journeys a pleasure to ride. The Cityliner scooters boast full lighting and indicators, as well as rear braking lights for greater safety whilst travelling in poor conditions. The Cityliner 312 and 412 combine excellent safety with performance, with new product features giving them the edge and advantage in this class.

MEYRA
ORTOPEDIA
We move people.

MODEL NO. 2.363/2.364

Cityliner₃₁₂
Cityliner₄₁₂

- 12 km/h (7,5 mph)
- Range up to max. 35 km (22 miles)
- User weight up to 150 kg (23,5 stone)
- Climbing ability up to 27% (with 75 kg/12 stone user weight)
- Motor power 400 W
- Smooth running, powerful motor
- Captain seat with adjustable headrest
- Illuminated, multifunctional LCD-display
- All round suspension
- 6-A-Charger
- Shopping basket
- Photosensor activates the lights when darkness increases
- Aluminium wheel rims



Electrical seat lift



Dimensions/size for transport

Length	120 cm/122 cm * (47,2/48" *)
Width	64 cm (25,2")
Seat height, adjustable	45 – 55 cm (17,7 – 21,7")
Seat depth	40 cm (15,7")
Seat width, adjustable	45 – 66 cm (17,7 – 25,9")
Height for transport (folded tiller, without seat)	65 cm (25,6")

Chassis

Automatic safety brake	•
Wheels	3/4 *
Wheel diameter front wheel	29 cm (11,4")
Wheel diameter rear wheel	29 cm (11,4")
Ground clearance	10 cm (3,9")
Obstacle height, max.	10 cm (3,9")
Turning radius	1,1 m/1,5 m* (43,3/59,1" *)

Electronic

Motor constant performance	400 W
Motor peak performance	1.600 W
Batteries	40 Ah (5 h) 50 Ah (20 h)
Controller	110 A

Performance

Travelling range, max. (under test conditions)	35 km (22 miles)
Technical climbing ability (with 75 kg/12 stone user's weight)	15° (27 %)
Safe ascending/descending gradient	10° (18 %)

Weights

Unladen weight	97,5 kg (106 kg) 215 lbs (233,7 lbs *)
Permitted user weight	130 kg (150 kg) **
Permitted total weight	20,5 stone (23,5 stone) ** 260 kg (41 stone)

Speed

Speed, max.	12 km/h (7,5 mph)
-------------	-------------------

Features / Options / Accessories

Colour	silver
Aluminium wheel rims	•
Lap belt	•
Lights	•
Brake light	•
Shopping basket, additional	•
Electrical seat lift **	•
Cane holder	•
Integrated bumper	•
Rear-view mirror, pair	•
Seat sliding rails	•

• Standard • Option

* Cityliner 412

** Safe ascending/descending gradient max. 12% (thus does not comply with EN 12184) and lower performance values

cbb 03.07 Technical product data subject to change. Colours may vary from those shown in brochures. Printing errors excepted. ID no. 261 510 901