

Cityliner<sub>315</sub>

Cityliner<sub>415</sub>



## CITYLINER 315/415 EXCLUSIVELY IN THE LEAD

The Cityliner 315 and 415 models are characterised by an integral bumper and ergonomically shaped steering column. Not only the design is exclusive, the equipment is special as well: cruise control, illuminated console, storage space, mirror and multi-function LCD display. Features such as the handbrake designed as a dual drum brake, lap belt and brake light combined with their rugged construction make these scooters your reliable companions. Guaranteed to make travelling a pleasure.

**MEYRA**  
**ORTOPEDIA**  
We move people.

# MODEL NO. 2.663/2.664

Cityliner<sub>315</sub>  
Cityliner<sub>415</sub>

- 15 km/h (10 mph)
- Range up to max. 50 km (31 miles)
- User's weight up to 200 kg (31,5 stone)
- Climbing ability up to 27% (with 75 kg/12 stone user's weight)
- Motor power 650 W
- LCD display shows time, date, etc.
- Luxury seat
- Speedometer/mileage indicator
- With shopping basket
- Lights comply with most road traffic regulations incl. indicators and hazard warning flasher
- Aluminium rims
- Suspension



Also as a three-wheel version, available from April 2006

## Dimensions/size for transport

Length	138 cm (54,3")
Width	68 cm (26,8")
Seat height, adjustable	44 – 54 cm (17,3"–21,3")
Seat depth	50 cm (19,7")
Seat width	50 cm (19,7")
Height for transport (steering column folded down, seat removed)	70 cm (27,6")

## Running gear

Automatic safety brake	●
Wheels	4
Steering wheel diameter	31,5 cm (12,4")
Driving wheel diameter	31,5 cm (12,4")
Ground clearance	11 cm (4,3")
Obstacle clearance, max.	10 cm (3,9")
Turning circle (radius)	1,7 m (66,9")

## Electronics

Motor continuous power output	650 W
Motor peak power	2.500 W
Batteries	64 Ah (5 h) 75 Ah (20 h)
Electronics	130 A

## Performance

Travelling range, max. (under test conditions)	50 km (31 miles)
Technical climbing ability (with 75 kg/12 stone user's weight)	15° (27 %)
Safe ascending/descending gradient (generally permissible)	10° (18 %)

## Weight

Unladen weight	163,5 kg (359,7 lbs)
Permissible user's weight	180 kg (200 kg)**
	28 stone (31,5 stone) **
Permissible total weight	375 kg (59 stone)

## Speed

Speed, max.	15 km/h (10 mph)
-------------	------------------

## Equipment features/options

Colour	silver
Aluminium rims	●
Lights	●
Brake light	●
Shopping basket, additional	●
Throttle pedal	●
Holder for mobility aids	●
Speed reduction on bends	●
Luxury seat	●
Fender	●
Rear-view mirror, left	●
Seat adjustment bar	●
Speedometer/mileage indicator	●
Warning panel	●

● Standard ● Optional

\*\* Safe ascending/descending gradient max. 15% (thus does not comply with EN 12184) and lower performance values

A.M. Mobility  
The Solihull Mobility Centre  
767-769 Old Lode Lane  
Solihull  
West Midlands  
B92 8JE

**MEYRA**  
**ORTOPEDIA**  
We move people.