



ORTOCAR
315 SP

ORTOCAR
415 SP



ORTOCAR 315 SP/415 SP NOTHING BUT THE BEST WITH THESE TOP MODELS

The Ortocar 315SP and 415SP meet all requirements. Both scooters have the "made in Germany" mark because the Meyra/Ortopedia Group manufacture them both. That means they cater for individual requirements and bespoke ordering is possible. With a quality that demands to be seen and experience in the build, only the best materials and evenly coloured plastics are used. New features include; larger wheels, front drum brakes, integral mudguards, powerful 80 Ah batteries and optimised standard equipment. Only the best for frequent travellers!!

MEYRA
ORTOPEDIA
We move people.

MODEL NO. 2.563/2.564

ORTO | CAR
315 SP

ORTO | CAR
415 SP

- 15 km/h (10 mph)
- Range up to max. 55 km (35 miles)
- User's weight up to 150 kg (23,5 stone)
- Climbing ability up to 30% (with 75 kg/12 stone user's weight)
- Motor power 700 W
- Adjustable sprung chassis
- Luxury seat
- Speed reduction on bends
- Slope warning
- With shopping basket
- Lights comply with most road traffic regulations including indicators and hazard warning flasher
- Speedometer/mileage indicator
- Aluminium rims
- Suspension



Dimensions/size for transport

Length	136 cm/143 cm *(53,5"/56,3"*)
Width	63 cm (24,8")
Seat height, adjustable	44 – 54 cm (17,3" – 21,3")
Seat depth	45 cm (17,7")
Seat width	50 cm (19,7")
Height for transport (steering column folded down, seat removed)	70 cm (27,6")

Running gear

Automatic safety brake	●
Wheels	3/4*
Steering wheel diameter	31,5 cm/36 cm (12,4"/14,2")
Driving wheel diameter	36 cm (14,2")
Ground clearance	13 cm (5,1")
Obstacle clearance, max.	10 cm (3,9")
Turning circle (radius)	1,6 m/1,7 m * (62,9"/66,9"*)

Electronics

Motor continuous power output	700 W
Motor peak power	2.800 W
Batteries	60 Ah (5 h) 80 Ah (20 h)
Electronics	130 A

Performance

Travelling range, max. (under test conditions)	55 km (31 miles)
Technical climbing ability (with 75 kg/12 stone user's weight)	17° (30 %)
Safe ascending/descending gradient (generally permissible)	11° (20 %)

Weight

Unladen weight	142/148 kg (312,4 / 325,6 lbs)*
Permissible user's weight	150 kg (23,5 stone)
Permissible total weight	310 kg (49 stone)

Speed

Speed, max.	15 km/h (10 mph)
-------------	------------------

Equipment features/options

Colour	silver
Aluminium rims	●
Amputee leg supports	●
Trailer	●
Trailer towing coupling	●
Lap belt	●
Leg supports (foot rest)	●
Lights	●
Indicators	●
Brake light	●
One-hand operation	●
Shopping basket, additional	●
Throttle pedal	●
Twist-grip throttle	●
Holder for mobility aids	●
Speed reduction on bends	●
Handle extension	●
Horn, amplified	●
Luxury seat	●
Fender	●
Walking frame holder	●
Rear-view mirror, left	●
Rear-view mirror, right	●
Seat adjustment bar	●
Slope warning	●
Speedometer/mileage indicator	●
Fixing device for transport	●
Warning panel	●

Standard
A.M. Mobility
 The Solihull Mobility Centre
 767-769 Old Lode Lane
 Solihull
 West Midlands
 B92 8JE

MEYRA
ORTOPEDIA
 We move people.